



"JOYA" COLLYRIUM APPLICATOR

GENERAL DESCRIPTION OF THE PRODUCT:

Device for making easy to apply drops in a patient eyes.
Its anatomical design, allows it for being holded with index and thumb fingers.
It has a central hole and two lateral holes in each side.

Composition

PP (Polypropylene)

Presentation

Plastic bag + cardboard box + instructions leaflet (1 u.)

INDICATIONS:

Suitable to apply collyrium oneself to its own eye without using any mirror, in so doing, we can decrease injuries risk produced by the contact of any collyrium bottle with ocular surface.

INSTRUCTIONS FOR USE:

You can see step by step the correct use of the collyrium applicator on the instructions leaflet.
Through images on the instructions leaflet, your can see its correct use according the different existent collyrium bottles.

MAINTENANCE:

Wash it after its use.
Do not share it with other persons.
Do not use it into oven, microwaves and do not boil it.

ROSER HUARTE PRATS
Technical Director
Biologist bar assn. No. 14413-C